


Application lab report

Preparation/ electrolytic polishing of aluminum samples



QATM-Preparation method

Objective:


This report shows the electrolytic polishing of unmounted aluminum samples (Figure 1) with our device **Qetch 1000**. The weld seam area should be made visible in order to be able to evaluate the weld seam. The preparation steps in gray fond are only a recommendation for cutting, mounting and grinding/polishing.




Figure 1: As-received samples.

QATM-Preparation method

Cutting










 Cutting			
Device	Cut-off disc	Anti-corrosion coolant	Clamping tool
Qcut 250 A	NF-A-250SC, 95012531	ATM CoolCut, 95004146	- Qtool 80 Vario Z2231200
Cutting method			
Automatic horizontal cutting method (with X-Axis)			
Parameters			
Feed speed	Pulse parameter	Cut-off disc rotational speed	
0,3 mm/s	Without pulse	3015 RPM	
Notes			

Cold mounting

 Mounting				
Consumable	Mixing ratio Volume	Curing time	Mold	Additional equip.
KEM 20	Powder : Liquide 2 : 1	15 Minutes	PP, ø 40 mm	- Dosing spoon - Mixing cup - Mixing stick
Notes				
- we recommend to mount the sample/s in the way that the sample backside is not covered with the mounting resin. The free backside of the sample will allowed direct contact of the mounted aluminum piece				

QATM-Preparation method

Grinding/Polishing

Device	Samples holder	Pressure mode				
Qpol 300 A1 Eco+	Z5446023	Single				
Step	MEDIUM		 RPM	 N	 min	
 Planar grinding	SiC-Paper P320	H ₂ O	150	100 ◀ ▶	15	Till flat (0:30 – 0:45)
 Grinding	SiC-Paper P600	H ₂ O	150	100 ◀ ▶	15	0:45
 Grinding	SiC-Paper P1200	H ₂ O	150	100 ◀ ▶	15	0:45
 Polishing	SIGMA	Diamond susp. Wb. Poly, 3 µm + wb. Lubricant	150	100 ▶▶	30	6:00
 Fine polishing	OMEGA	Eposil F 0,1 µm	100	90 ◀ ▶	20	4:00 (Rinsing time: 10 s)

Notes

- Pre-dosing for 3 µm: 3 s
Dosing interval and dosing duration for diamond susp. 3µm:
Every 35 s for 1,5 s
- Dosing interval and dosing duration for lubricant:
Every 75 s for 1,5 s
- Dosing interval and Dosing duration for fine polishing suspension:
Every 20 s for 1,0 s

QATM-Preparation method


Electrolytic polishing and etching

Device					
--------	--	--	--	--	--

Qetch 1000

Step	Electrolyte	Voltage [V]	Pumping power [%]	Mask [cm ²]	Duration [min]
------	-------------	-------------	-------------------	-------------------------	----------------

Electrolytic polishing

	Polishing (electrolytic)	Electrolyte K1	10	50	1,5	0:20
---	--------------------------	----------------	----	----	-----	------

Electrolytic Etching

	Etching (electrolytic)	-	-	-	-	-
---	------------------------	---	---	---	---	---

Notes - before using the Qetch 1000 the sample/s must be pre-grind with SiC-Paper, P600

QATM-Preparation method




Figure 2: Overview of the sample placed on the mask.




Figure 3: Sample after the electrolytical polishing process with visible indent area.

QATM-Preparation method



Figure 4: Sample A – like we have received – 25 x




Figure 5: Sample A - after Qetch 1000 – 10V 20s (pre-grinded with SIC 600) – 25 x

QATM-Preparation method




Figure 6: Overview of sample B after electrolytical polishing process with Qetch 1000




Figure 7: Sample B - after Qetch 1000 - 10V 20s (pre-grinded with SIC 600) - 25 x

QATM-Preparation method




Figure 8: Sample C – like we have received – 25 x




Figure 9: Sample C - after Qetch 1000 – 10V 20s (pre-grinded with SIC 600) – 25 x

QATM-Preparation method




Figure 10: Sample D – like we have received – 25 x




Figure 11: Sample D - after Qetch 1000 – 10V 20s (pre-grinded with SIC 600) – 25 x

QATM-Preparation method

Conclusion:

The Qetch 1000 can be used very well for the preparation of the weld seam area. The weld seam is very clearly visible. A macroscopic analysis could be performed. Before starting the process with the Qetch 1000 – the sample must be pre-grind with SIC-Paper, P600.

The preparation steps marked in grey are recommendations without having carried them out.