



电解抛光腐蚀机

QETCH 100 M

Qetch 100 M是一款便携式由电池供电的电解抛光腐蚀机。它结构紧凑坚固，是现场使用的理想选择。

其操作简便性和功能强大性独一无二。带有树脂涂层的铝制外壳满足功能性、稳定性和设计的最高要求。所有敏感部件出于运输和现场操作的考虑都得到很好的保护。为了最大程度的舒适，该机还配备了手柄和背带。

优点

- | 便携易用的电解抛光和腐蚀单元
- | 微处理器控制的全自动程序运行
- | 使用锂电池即可运行机器，无需电源
- | 集成抛光单元，开始/停止功能
- | 彩色LED显示屏显示当前的工作状态
- | LCD显示屏显示了所有重要的工作状态
- | 电压，电流和腐蚀/抛光时间可调

电解抛光腐蚀机 QETCH 100 M

技术参数

最大功率	160 W
最大电压	90 V DC
最大流量	2A
抛光面积	~ \varnothing 9 mm
宽 x 高 x 深	380 x 355 x 180 mm
电解液槽	1升
重量	~ 8.5 kg (不含电解部分)

www.qatm.cn/qetch-100-m

ORDER DATA

BASIC EQUIPMENT

MACHINE (ALL ITEMS REQUIRED)

M6201000  Qetch 100 M (Kristall 650) Base machine

POWER SUPPLY (1 ITEM REQUIRED)

A6201000 Qetch 100 M power supply 220-230V/50/60Hz

A6201001 Qetch 100 M power supply 110V/60Hz (1Ph/N/PE)


ELECTRICAL BOARD FOR PORTABLE GRINDER (TO BE ORDERED WITH THE MACHINE)

A6201004 Unit for using without portable grinder

A6201005 Unit prepared for using portable grinder

PORTABLE GRINDER UNIT FOR QETCH 100 M

Z6201001  Accessory-Set Qetch 100 M

95003955  Plastic Bottle 1000ml

Z6201000 Handle holder